

Heron has partnered with major full-line wholesalers and specialty distributors for simplified, fast, and cost-effective ordering. To order ZYNRELEF, contact your preferred distributor. If you have any distribution-related questions, please contact your account manager or call Heron Connect at 1-844-HERON11 (1-844-437-6611).

DISTRIBUTOR NAME	WEBSITE	PHONE
Hospitals		
AmerisourceBergen Drug Corporation	www.amerisourcebergen.com	1-844-222-2273
Cardinal Health, Inc	www.cardinalhealth.com	1-800-926-3161
McKesson Pharmaceutical	connect.mckesson.com	1-855-625-6285
Morris & Dickson Co LLC	www.morrisdickson.com	1-800-388-3833
Ambulatory Surgery Centers (ASCs)		
Besse Medical	www.besse.com/home	1-800-543-2111
Cardinal Health SPD	www.cardinalhealth.com/spd	1-866-677-4844
CuraScript SD	www.curascripts.com	1-877-599-7748
Henry Schein	www.henryschein.com/medical	1-800-772-4346
McKesson Medical-Surgical	mms.mckesson.com	1-855-625-4677
Metro Medical	www.metromedicalorder.com	1-800-768-2002

Package Information, Storage, and Handling

ZYNRELEF is supplied as a kit consisting of a nonsterile single-dose glass vial (containing sterile active ingredients) and sterile components needed for preparation and application. Store ZYNRELEF kits at controlled room temperature (20°C to 25°C [68°F to 77°F]). Protect from moisture and light.

Product Code (NDC)	Bupivacaine/Meloxicam	Volume	Kit Contents (in addition to single-dose vial)
47426-301-02	400 mg/12 mg	14 mL	1 vented vial spike, 2 x 12-mL Luer lock syringes, 2 Luer lock cone-shaped applicators, 2 syringe tip caps
47426-303-01	200 mg/6 mg	7 mL	1 vented vial spike, 1 x 12-mL Luer lock syringe, 1 Luer lock cone-shaped applicator, 1 syringe tip cap
COMING SOON 47426-501-02	400 mg/12 mg	14 mL	1 x 20-mL sterile Vial Access Needle (VAN), 2 x 12-mL Luer lock syringes, 2 Luer lock applicators, 2 syringe tip caps
COMING SOON 47426-503-01	200 mg/6 mg	7 mL	1 x 10-mL sterile Vial Access Needle (VAN), 1 x 12-mL Luer lock syringe, 1 Luer lock applicator, 1 syringe tip cap

Note: VAN kit dimensions: 6.25" (length) x 2.00" (height) x 4.00" (width). All other kit dimensions: 4.04" (length) x 1.60" (height) x 7.04" (width).

HeronConnect Provides Extensive Access and Support

Replacement parts

Prepare and administer ZYNRELEF using only parts provided in the ZYNRELEF kit. If parts are missing or damaged, Heron will provide a replacement kit or component, including:

- **Luer lock cone-shaped applicators** (GTIN 00347426902106)
- **Vented vial spikes** (GTIN 00347426903059)
- **12-mL Luer lock syringes** (GTIN 00347426904087)

Billing/coding support

Dedicated Reimbursement Counselors will ensure you have the information needed to bill and code ZYNRELEF. Contact Heron Connect at **1-844-HERON11 (1-844-6611)** 8 AM to 5 PM ET, Monday through Friday. For more information, visit [HeronConnect.com](https://www.HeronConnect.com)



Please see accompanying full Prescribing Information, including Boxed Warning.